

# *ToxTrac*: a fast and robust software for tracking organisms

## User guide for version 2.61

(Submitted to the Journal Methods in Ecology and Evolution)

Alvaro Rodriguez<sup>1</sup>, Hanqing Zhang<sup>1</sup>, Jonatan Klaminder<sup>2</sup>, Tomas Brodin<sup>2</sup>, Patrik Andersson<sup>3</sup>, Magnus Andersson<sup>1</sup>

<sup>1</sup>Department of Physics, Umeå University, 901 87 Umeå, Sweden, <sup>2</sup>Department of Ecology and Environmental Science, Umeå University, 901 87 Umeå, Sweden, <sup>3</sup>Department of Chemistry, Umeå University, 901 87 Umeå, Sweden

**Keywords:** animal behavior; guppy; salmon; zebrafish; tadpole; cockroach; Kalman filter; ecology; tracking software

## Contents

Starting Guide.....	5
1. Requirements.....	5
2. Installing and running the software .....	5
3. Known Issues .....	5
4. Recording .....	6
Video Files .....	6
Tracking areas and background.....	6
Animals .....	6
5. Experimental Setup .....	6
Arena.....	7
Illumination .....	7
Basic Functions .....	9
1. Start Screen.....	9
2. Project Files.....	9
Project Overview.....	11
1. Project overview screen.....	11
2. Calibration algorithm .....	12
Camera model.....	12
Calibration using patterns.....	12
Manual calibration.....	13
3. Calibration screen.....	14
Arena definition.....	15
1. Arena definition algorithm .....	15
Automatic selection .....	15
Manual selection.....	16
2. Arena definition screen.....	17
Automatic selection .....	18

Manual selection.....	19
Detection .....	20
1. Detection algorithm.....	20
Background subtraction.....	20
Animal detection .....	21
Filtering .....	21
2. Detection screen .....	22
Tracking .....	23
1. Kalman tracking algorithm.....	23
Feature Extraction .....	24
Acceptance conditions.....	25
Collision .....	26
Delete conditions.....	26
2. Fragment Identification .....	26
3. Tracking screen .....	27
Results .....	28
1. Statistics screen .....	28
Individual statistics.....	28
Population statistics .....	29
2. Graphical outputs screen .....	29
Individual graphical outputs .....	29
Population graphical outputs .....	31
3. Output Files .....	31
Tracking ( <i>Traking.txt</i> ) .....	33
Tracking in Real Space coordinates ( <i>Tracking_RealSpace.txt, Trajectory.jpeg</i> ).....	33
Instantaneous Speed ( <i>Instant_Speed.txt</i> ) .....	35
Instantaneous Acceleration ( <i>Instant_Accel.txt</i> ) .....	35

Distance to Edges ( <i>Dist_Edges.txt</i> , <i>Dist_Edge_[edge].jpeg</i> ).....	36
Exploration ( <i>Exploration.txt</i> , <i>Exploration.jpeg</i> ).....	38
Distance to Mean Position ( <i>Dist_MeanPos.txt</i> , <i>Dist_Mean_Pos.jpeg</i> ).....	41
Distance to Center Position ( <i>Dist_CenterPos.txt</i> , <i>Dist_Center_Pos.jpeg</i> ).....	43
Transitions ( <i>Transitions.txt</i> ) .....	45
Frozen Events ( <i>FrozenEvents.txt</i> ) .....	45
Stats ( <i>Stats.txt</i> ).....	45
Advanced parameters .....	46
1.    Graphic parameters ( <i>ColorIni.txt</i> ) .....	46
2.    Configuration parameters ( <i>Configuration.txt</i> ).....	47
References .....	54

## Starting Guide

### 1. Requirements

*ToxTrac* has been developed for Windows in C++ using Visual Studio 2015 with the add-on Qt5Package. In addition we have used OpenCv3.0, an open-source computer vision library available at <http://opencv.org>, Qt5.6.0 open-source, a library for building interfaces available at <https://www.qt.io>, and LibXL, a library for writing and reading excel files available at <http://www.libxl.com>.

*ToxTrac* requires Windows 7 or later, and is created for 64-bit hardware. We recommend a minimum of 8 GB of RAM memory and enough hard drive free space to handle all recorded video files. A 2.0+ GHz Quad core or higher is recommended for proper performance.

### 2. Installing and running the software

The *ToxTrac* project, containing the latest software version, the documentation of the program and other resources is hosted at <https://toxtrac.sourceforge.io>.

To install *ToxTrac*, it is only necessary to execute the .exe Windows installer file provided.

Visual Studio 2015 or Visual C++ x64 2015 Redistributable Packages (or newer) are required to run *ToxTrac*. The installer will automatically download and install this component if necessary.

Once Visual Studio or the Visual C++ Redistributable Package is installed, follow the instructions of the install process. *ToxTrac* and all necessary components will be installed in the selected folder and a shortcut in the Desktop and the start menu will be created. Now *ToxTrac* is ready.

An updated codec pack is recommended to properly handle video files with *ToxTrac*. We have used K-Lite Codec Pack, available at <https://www.codecguide.com>.

It is also recommendable to install any spreadsheet package compatible with .xls (Microsoft Excel files) to view the statistical file generated by *ToxTrac*.

To inspect and edit video files we recommend VirtualDub, available at <http://virtualdub.org>.

To inspect and edit image files, we recommend GIMP, available at <https://www.gimp.org>.

### 3. Known Issues

- A minimum screen resolution of 1280x960 is required to display properly the interface windows.
- *ToxTrac* is aimed to work with 96 dots per inch screens (the Microsoft Windows operating system default display). Changing the dpi, for example to increase the font size in the windows accessibility configuration, will not increase the text size. Instead, *ToxTrac* will try to rescale all text elements to avoid clipping. This behavior may result in difficulties to visualize the interface for some screens, especially in tablets or some portable devices.
- *ToxTrac* does not recognize non-standard folder of file names. Though spaces are permitted, the program will now work properly if: The used file names, the project, the calibration or the video folders contain non English characters.
- Using the fragment identification algorithm may require a lot of RAM memory. If the used memory exceeds the system capacity, *ToxTrac* may crash. The fragment identification algorithm can be configured to avoid this issue, (for example using a TCM depth of 8 pixels, reducing the TCM radius or the maximum number of samples per fragment). Additionally, changing the maximum number of fragments will force the program to free the memory when the number of trajectory fragments in memory reaches a certain limit. The identification algorithm can be partially or totally deactivated to reduce computational time and prevent any memory problem.

- If the project location is in a folder without write permission, *ToxTrac* will not be able to save the project data and will show a warning video.

## 4. Recording

### Video Files

*ToxTrac* supports a wide variety of .avi video files with any resolution and framerate, including MPEG-4 and x264 compressed .avi video files. However, the video should of course have the highest possible quality.

If the video is cut in several files, all the pieces should be placed in the same folder and file names should end in a correlative number according to the file order (example: video01.avi, video02.avi, ...).

Multiple video sequences (each sequence composed by multiple video files) can be analyzed at the same time (if experimental conditions are not changed).

Video resolution should be high enough so that the animal size is at least 50 pixels, and framerate should be high enough so the animal area in consecutive frames overlaps. We find that 25fps is enough for most experiments, but with fast moving animals, and especially in multiple animals experiment, a higher framerate may be advisable.

*ToxTrac* will automatically convert the video image format to a grayscale image, so color information is not needed.

### Tracking areas and background

*ToxTrac* will detect and track animals in rectangular pieces of the image containing the arenas where we want to observe the animals. Inside the arena, the tracking areas are defined as uniform bright regions with no particular shape, where the tracked objects can be detected. If the arenas have dark corners or edges, these should be excluded from the tracking area.

Ideally, the background color in the tracking areas should as homogeneous as possible, and brighter than the animal, with the highest possible contrast.

The presence of different background objects (dark objects appearing inside the tracking areas) is acceptable, and will be managed by the system in different ways: Static or moving objects much smaller than the animals will not affect segmentation, because they can be filtered out. Static objects of any size can be removed using the background subtraction technique (this may cause animals not to be detected until they start to move). Alternatively, the arena selection tool can be used to exclude parts of the background.

Objects which cannot be separated from the animals (moving objects, or static objects which are not excluded from the arena) will difficult or impede the tracking process. It is especially important to use a background as free of reflections or shadows as possible (see lightning) and to exclude the dark edges of the experimental setup from the background (see arena definition).

### Animals

Multiple animals can be tracked in the same arena, however, we don't recommend using more than 10-20 animals in a single experiment if it is important to keep the identity of the animals during the experiment. This value however depends on the occlusion degree observed in the video.

## 5. Experimental Setup

The experimental setup is typically formed by the arenas, the camera and the illumination elements (lights, filters, diffusers...). During an experiment, the experimental setup should be isolated from external interference and external light variations.




Figure 1. Arena schematic.

It is very important that the experimental setup present the same conditions and it is not moved during an experiment, between the calibration and the experiment, or between experiments, if the same calibration data is used or if they are going to be processed together.

### Arena

The walls of the arena should not cast strong shadows or reflections. This can be achieved by studying the camera position, the lightning conditions and the arena materials. Transparent or translucent walls will not cast shadows, but may have reflections; and opaque walls will not have reflections but may cast shadows. Also wall height can be adjusted to reduce these effects, and walls can be in most cases excluded from the tracking area.

### Illumination

Without an appropriate illumination, the task of tracking is impossible. The light conditions determine the type, position, angle and intensity of the beams incident in the arenas which will be then registered by the camera.

We recommend two types of illumination to be used in the tracking system:

- Diffuse illumination, this light preserves the texture details and mitigates shadows.
- Backlight illumination, this light will highlight the shapes and will not cause shadows. However, it will also hide the textures of the objects.

A scheme of how to construct these illumination types is shown in Figure 2, and Table 1 shows how each type of light interacts with different characteristics of the objects.




Figure 2. Dark filed and backlight illumination.

Table 1. Dark filed and backlight illumination characteristics.

Feature	Backlight	Dark Field
Absorption (Changes in light absorption from the object)	None	Minimal effect
Texture (Changes un surface texture)	None	Textured surfaces brighter than polished
Elevation (Changes in height from surface to camera z axis)	None	Outer edges are bright
Shape (Change on shape or contour along x/y axis)	Shows outside contours	Contours highlighted, flat surfaces darker than raised
Translucency (Changes in density-related light transmission)	Shows changes in translucency vs. opaqueness	None

## Basic Functions

### 1. Start Screen




Figure 3. Start Screen.

- 1: Access to the load/save menu, allowing different options: “**New Project**” (deletes all current data, and reloads color and configuration files), “**Save Project**”, “**Load Project**” or “**Merge Project**” in the current one (merges another project in the current one, this requires that both projects are recorded in the same conditions, with the same calibration data and the same arena definition parameters).
- 2: Settings menu, this allows the user to navigate through the different settings of the software. These are divided into four different panels: “**Project Overview**”, “**Arena Definition**”, “**Detection**”, and “**Tracking**”.
- 3: Results menu, (only available after analyzing the video sequences), this allow the user to see the statistical data of the video sequences and the graphical outputs.
- 4: Tracking/Stop, Process the videos with the current configuration. It automatically saves the project when the analysis starts and when it finishes. During the analysis, this button will allow the user to stop the current processing, in a safe point in the body tracking or in the fragment identification algorithm. If a process is stopped, all progress is lost.
- 5: Window Menus, Standard windows buttons to minimize, maximize and close *ToxTrac* main window.

### 2. Project Files

A project is a coherent tracking analysis, identified by the project name (*out.pnam*). One unique calibration and configuration will be used in one project. One project can have several video sequences and one video sequence can have several video files. All sequences are expected to be recorded in the same experimental conditions and with the same camera parameters.

When pressing the Analyze button, every sequence in the project will be analyzed, and the results will be located in a set of subfolders in the project folder. Also, all results will be joined to estimate the statistics of the entire population analyzed.

Table 2. Project files

Project File	Example of the File	Description
pname.tox	C:/ProjectFolder/pname/pname_Input.txt C:/ProjectFolder/pname/pname_Configuration.txt C:/ProjectFolder/pname/pname_Arena.txt C:/ProjectFolder/pname/pname_ArenaNames.txt C:/ProjectFolder/pname/pname_Calibrator.txt C:/ProjectFolder/pname/pname_Output.txt	Main project file. Contains a link to Input, output and configuration files Is generated by the application when saving or analyzing, and can be loaded in the application Can be added to another compatible project, joining both populations.
pname_Input.txt	1 2 C:/VideoFolder/Seq1_0.avi C:/VideoFolder/Seq1_1.avi 0 2054	Video input file. Contains the information of the video sequences. For each video sequence: the corresponding video files A video and frame number is stored to define a reference frame where arenas will be defined.
pname_Configuration.txt	CALIBRATION_PARAMETERS [...] ARENA_DEFINITION_PARAMETERS [...] BACKGROUND_PARAMETERS [...] PREPROCESSING_PARAMETERS [...] DETECTION_PARAMETERS [...] KALMAN_FILTER_PARAMETERS [...] KALMAN_MULTITRACKING_PARAMETERS [...] DATA_ANALYSIS_PARAMETERS [...] MAIN_VIDEO_PARAMETERS [...]	Configuration file. Contains all the main parameters used in the program. Some of these parameters are accessible in the interface. These parameters can be modified in the text file, to modify the algorithms behavior.
pname_Arena.txt	1 200 4 1522 1075	Arena definition file Contains the corner coordinates in pixels, of the manually drawn arenas.
pname_ArenaNames.txt	1 Arena1	Arena Names file Contains the names of the arenas introduced by the user
pname_Calibrator.txt	[...] 1 1 0 0 10 0 0 0 0 9.5 0 1 0 0 1 [...]	Calibration parameters file. Contains the camera matrix, the rotation matrix, the displacement vectors and the distortion coefficients of the calibration. This file is used to scale the primary output and also to remove the distortion from every video frame.
pname_Output.txt	1 2 C:/ProjectFolder/ Seq1 / Tracking_0.txt C:/ProjectFolder/ Seq1 / Tracking_1.txt	Primary output file. Points to the tracking results. All results will be extracted from these files by the application. This output will generate a file for each arena and sequence.

## Project Overview

### 1. Project overview screen

This screen shows the selected videos to track, and allows to add, remove, or reorder them. In addition, the user can modify basic parameters.




Figure 4. Project overview screen.

- 1:** Load a new video sequence in the list. To add a sequence formed of multiple videos named with correlative numbers, select the first video of the sequence and all files representing video fragments of the same video will be added to the list as a video sequence. Each sequence will be analyzed as one entity.
- 2:** Select project location in the hard drive, all results will be saved in this folder, be sure you have writing permissions in this folder.
- 3:** Restart the project, removes all results, videos, and calibration data, keeping current configuration.
- 4:** Remove all results computed for current videos.
- 5:** Open calibration window.
- 6:** Select the starting point in mins (for all sequences) for the analysis.
- 7:** Select the ending point in mins (for all sequences) for to analysis.
- 8:** List of videos, right click in one of the videos allows to move up or down, or to delete the corresponding sequence in the list. It also shows the status of the current video, and the completion rate of the different stages of the process.
- 9:** Interpolate holes in the trajectory, using a linear interpolation algorithm. This parameter can be changed without redoing the analysis.
- 10:** Maximum size of a trajectory hole (in frames) where interpolation will be applied. This parameter can be changed without redoing the analysis.
- 11:** Use a moving average to smooth trajectories. This parameter can be changed without redoing the analysis.

## 2. Calibration algorithm

### Camera model

In practice, due to small imperfections in the lens and other factors, some distortions are integrated into the image. These distortions can be modeled using the following parametric equations<sup>2</sup>:

$$d(x) = x \frac{1 + k_1 r^2 + k_2 r^4 + k_3 r^6}{1 + k_4 r^2 + k_5 r^4 + k_6 r^6} + 2 p_1 xy + p_2 (r^2 + 2x^2) + s_1 r^2 + s_2 r^4, \quad (1)$$

$$d(y) = y \frac{1 + k_1 r^2 + k_2 r^4 + k_3 r^6}{1 + k_4 r^2 + k_5 r^4 + k_6 r^6} + p_1 (r^2 + 2y^2) + 2 p_2 xy + s_3 r^2 + s_4 r^4, \quad (2)$$

$$r^2 = x^2 + y^2, \quad (3)$$

where  $x$  and  $y$  are spatial coordinates,  $r$  is the distance to the lens optical center,  $d(x)$  and  $d(y)$  are the corresponding distorted coordinates,  $k_i$  are the radial distortion coefficients,  $p_i$  are the tangential distortion coefficients and  $s_i$  the prism distortion coefficients.

The distortion coefficients do not depend on the scene viewed, thus they also belong to the intrinsic camera parameters. And they remain the same regardless of the captured image resolution.

### Calibration using patterns

Calibration is performed using a sequence of images of a calibration pattern. The calibration images should be recorded using the same conditions as in the intended experiment and in the same image resolution. The calibration pattern has a shape of a black and white chess-board pattern, without edge lines and with white edges. It must be printed in high resolution, and can be automatically generated in some image editing tools such as GIMP. The images should be taken in a well illuminated environment, in the plane of interest of the experiment, and using the same camera position and camera parameters as in the actual experiment. The calibration images should be placed in the same directory and with a file name ending in correlative numbers. An example of calibration pictures can be seen as follows.




Figure 5. Example of calibration views taken for a 4 arena experiment.

The calibration proceeds as follows:

The user first selects the first image of the sequence with *ToxTrac*.

The user selects one of the calibration images to define the camera pose.

A corner detection algorithm finds the positions of all the squares in the images. If the found positions do not correspond with the number of rows (*cal.rows*) or columns (*cal.cols*) of the pattern for a particular image, the algorithm will return an error, and that image should be eliminated from the sequence.

The algorithm will use the known dimensions of the squares of the calibration pattern (*cal.size*) and their positions in the image to estimate the calibration parameters using the global Levenberg-Marquardt optimization algorithm.

The user can select different distortion models (*cal.dist*), using different subsets of the distortion coefficients from equations 8-10. Available models are: “**Radial 3**” (dist. coeffs.:  $k_1, k_2, k_3$ ), “**Radial 3 + Tangential 2**” (dist. coeffs.:  $k_1, k_2, k_3$ ), “**Radial 6 + Tangential 2**” (dist. coeffs.:  $k_1, k_2, k_3, p_1, p_2, k_4, k_5, k_6$ ), “**Radial 6 + Tangential 2 + Prism 4**” (dist. coeffs.:  $k_1, k_2, k_3, p_1, p_2, k_4, k_5, k_6, s_1, s_2, s_3, s_4$ ).

### Manual calibration

When it is not possible to obtain views of the calibration patterns, calibration can be performed manually, introducing the parameters of a camera model.

The most straightforward way to do this is to assume that the alignment error and the distortion of the camera can be disregarded. In this case, translation vectors, distortion parameters, and the Euler pose rotation vectors should be set to 0. The calibration proceeds as follows:

- The user measures the horizontal and vertical mm to pixel scale of the image. For example a horizontal scale of 1:10 (1 mm = 10 pixel) and a vertical scale of 1:9.5 (1 mm = 9.5 pixel).

- The scale factors are introduced in the camera matrix parameters  $f_x$  and  $f_y$ . The camera matrix for the example values is shown in equation 11.

$$M = \begin{bmatrix} 10 & 0 & 0 \\ 0 & 9.5 & 0 \\ 0 & 0 & 1 \end{bmatrix}, \quad (4)$$

### 3. Calibration screen




Figure 6. Calibration screen, preview.

- 1:** Select a calibration sequence. The user must select the first the first image of the sequence. The calibration images must have a name ending with correlative numbers and should be placed in the same directory
- 2:** Select previous image for pose estimation.
- 3:** Select next image for pose estimation.
- 4:** Show/hide the advanced calibration parameters.
- 5:** Number of columns of the calibration pattern.
- 6:** Number of rows of the calibration pattern.
- 7:** Dimensions in mm of the squares of the calibration pattern.
- 8:** Distortion model.
- 9:** Try to estimate calibration parameters using the selected calibration file, and the current parameters.
- 10:** Reset the calibration parameters.
- 11:** Accept calibration.
- 12:** Cancel, closes the window, no changes will be made in the calibration data.
- 13:** View of the current image of the calibration sequence. If calibration if successful, it will show the detected features of the pattern.




Figure 7. Calibration screen, advanced.

- 1: Camera matrix, the pin-hole projection parameters of the camera.
- 2: The rotation matrix, is used to determine the orientation of the camera. Estimated for the current pose. Can be modified using the Euler angle representation (rotation angles) of the matrix.
- 3: Translation vectors used to determine the relative position of the camera Estimated for the current pose.
- 4: Distortion parameters, they model and correct the distortion caused by imperfections in the optical system of the camera.

## Area definition

An arena is a closed and controlled area, where a tracking experiment will take place, the arena is separated physically from the outside and from other tracking arenas.

Additionally, inside each arena, the tracking area must be defined prior to the execution. The tracking area should be a white uniform well illuminated area where the objects of interest will be studied. The objects of interest should appear as dark areas with high contrast inside the tracking area.

It is important than areas outside the white uniform region where the objects of interest are located are not inside the tracking areas or present also a white and uniform color. To void phenomena such as external objects, reflections, bubbles, shadows or interferences with walls and corners are ruled out.

The arena definition algorithm, provide a semiautomatic and easy to use tool to define and visualize the arenas and tracking areas.

### 1. Arena definition algorithm

#### Automatic selection

The arena selection algorithm takes a sample image of from an input video sequence to define the different tracking areas.

The algorithm to define arenas and tracking areas proceeds as follows:

- The user selects one frame of the video sequence to use as reference for the algorithm.
- Preprocessing: After the image has been obtained, the calibration model is used to create a distortion map for every pixel of the image, and then interpolation is used to create a distortion-free image. Then image is converted to a 8-bit grayscale, and normalized to value of 0-255.
- Segmentation: The objective of this operation is to separate the tracking areas, which by definition are uniform bright region of the images. First, an intensity value (*roi.thre*) is selected by the user to threshold the image into two binary sets. Finally, a closing mathematical morphological operation is executed. In mathematical morphology, the closing of a binary image  $A$  by a structuring element  $B$  is the erosion of the dilation of that set. This operation removes the holes and imperfections on the area selected by the previous operations. The size of the structuring element (*roi.elms*) and the iterations of the dilation (*roi.dilt*) and erosion (*roi.erot*) operations can be selected by the user.
- Arena and area creation: First, the areas obtained in the previous step are filtered according to its size (*roi.mins*). Then, the resulting areas which possess an arbitrary shape are approximated by a polygon so that the number of vertices is the smaller possible and the distance between vertices is less or equal to the precision specified by the user (*roi.poly*). This step simplifies the shape of the selected area, and it helps to eliminate some irregularities at the edges. These polygons will constitute the tracking areas. Finally, arenas are defined as the minimum rectangular regions in the image containing each one of the tracking areas.

In the application each arena will be defined as image region, containing a polynomial tracking area. Arenas are obtained by a map describing the undistorted positions of its pixels. Each arena constitutes unit processed in parallel with an independent tracking algorithm.

### Manual selection

The arena selection algorithm takes a sample image of from an input video sequence to define the different tracking areas.

The algorithm to define arenas and tracking areas proceeds as follows:

- The user selects one frame of the video sequence to use as reference for the algorithm.
- Preprocessing: After the image has been obtained, the calibration model is used to create a distortion map for every pixel of the image, and then interpolation is used to create a distortion-free image. Then image is converted to a 8-bit grayscale, and normalized to value of 0-255.
- Arena creation: The user manually draws an arbitrary number of rectangular shapes, which will constitute the arenas. Each arena constitutes unit processed in parallel with an independent tracking algorithm and each arena will contain a tracking area where the animals will be visible.
- Segmentation: First, an intensity value (*roi.thre*) is selected by the user to threshold the image into two binary sets. Finally, a closing mathematical morphological operation is executed. In mathematical morphology, this operation removes the holes and imperfections on selected area. The size of the structuring element (*roi.elms*) and the iterations of the dilation (*roi.dilt*) and erosion (*roi.erot*) operations can be selected by the user. The largest selected area inside each arena defined by the user is selected.
- Area creation: The largest selected area inside each arena defined by the user is selected. The user has the option to fit each of these areas to a circular shape (*roi.fite*), defined by the

minimum enclosing circle containing these areas, and the user can reduce the radius of the circle, by a selected number of pixels (*roi.redr*).

## 2. Arena definition screen




Figure 8. Arena definition screen.

- 1:** Select previous sequence. (Allows to preview the results of the current parameters in all the video sequences)
- 2:** Select next sequence.
- 3:** Select previous arena. The arena currently selected is displayed as a green rectangle, and other arenas are displayed in blue.
- 4:** Select next arena.
- 5:** Allows the user to introduce a name for the current arena.
- 6:** Select arena definition algorithm. (Automatic selection or manual selection).
- 7:** Open arena definition window, for the selected algorithm.
- 8:** View of the current arena selection. Shows a preview of the current video sequence and highlights the current arena. The tracking areas are displayed in red, and the arenas are shown as colored rectangles with the name over imposed.

## Automatic selection




Figure 9. Arena definition screen, automatic selection.

- 1: Threshold, this parameter selects the minimum intensity level (normalized to 1-255) of the tracking areas.
- 2: Polygon Flexibility, defines the maximum distance between vertices in the polygon with the smaller possible number of vertices used to define the tracking areas.
- 3: Select previous sequence. (Allows to preview the results of the current parameters in all the video sequences)
- 4: Select next sequence.
- 5: Select previous video of the current sequence.
- 6: Select next video of the current sequence.
- 7: Change the current frame of the selected video (when the button is pressed a new random frame of the video will be selected and displayed), the selected frame will be used for the selected sequence in the arena definition algorithm.
- 8: Minimum area, minimum number of pixels to constitute a tracking area.
- 9: Mask Size, size (diameter) of the structuring element in the closing mathematical morphology operation.
- 10: Dilation Iterations, number of dilation operation in the closing operation.
- 11: Erosion iterations, number of erosion operation in the closing operation.
- 12: View of the current arena selection. The tracking areas are displayed as red areas, and the arenas as blue rectangles.
- 13: Accept the arena selection.
- 14: Cancel, closes the window, no changes will be made.

## Manual selection




Figure 10. Arena definition screen, manual selection.

- 1:** Threshold, this parameter selects the minimum intensity level (normalized to 1-255) of the tracking areas.
- 2:** Select previous sequence. (Allows previewing the results of the current parameters in all the video sequences)
- 3:** Select next sequence.
- 4:** Select previous video of the current sequence.
- 5:** Select next video of the current sequence.
- 6:** Change the current frame of the selected video (when the button is pressed a new random frame of the video will be selected and displayed), the selected frame will be used for the selected sequence in the arena definition algorithm.
- 7:** Select previous arena. The arena currently selected is displayed as a green rectangle, and other arenas are displayed in blue.
- 8:** Select next arena.
- 9:** Remove selected arena.
- 10:** Add arena, opens a drawing window, which can be freely resized, and where the user can draw a rectangular area with the mouse. The drawing window can be closed pressing in the corner of the window or pressing enter.
- 11:** Remove all arenas.
- 12:** Mask Size, size (diameter) of the structuring element in the closing mathematical morphology operation.
- 13:** Dilation Iterations, number of dilation operation in the closing operation.
- 14:** Erosion iterations, number of erosion operation in the closing operation.
- 15:** Fits the tracking areas to their minimum enclosing circles.
- 16:** If the tracking areas are fitted to their minimum enclosing circles. The radius of the circles can be reduced by an arbitrary amount.

- 17:** View of the current arena selection. The tracking areas are displayed as red areas, and the arenas as blue or green rectangles.
- 18:** Accept the arena selection.
- 19:** Cancel, closes the window, no changes will be made.

## Detection

### 1. Detection algorithm

This algorithm will detect dark moving animals in a bright homogenous background. If after the arena selection procedure, the tracking area still contains static objects of a significant size, they can be removed using the background subtraction technique. However, this technique may only be used if the animals do not remain stationary (especially at the beginning of the video).

#### Background subtraction

The background subtraction process removes the static elements of the image, using a dynamic background modelling technique. According to this, every pixel of the scene must be matched to the background or foreground category. To this end a widely used model based in estimate the RGB color space probability distribution for every pixel in the image has been chosen. The segmentation algorithm works using Bayesian probability to calculate the likelihood of a pixel  $x_{ij}$ , at time  $t$  in coordinates  $(i,j)$ , being classified as background (*BG*) or foreground (*FG*). This is expressed as follows:

$$p(BG|x) = \frac{p(x|BG)p(BG)}{p(x|BG)p(BG) + p(x|FG)p(FG)}, \quad (5)$$

In a general case, we can assume that we don't know about the foreground objects and we may assume that  $p(BG)=p(FG)=0.5$  or we can use a different probability (*bgs.ratb*) according to the knowledge of the scene. The background model will be referred as  $p(x|BG)$  and we will decide that the pixel belongs to the background if  $p(x|BG)$  is higher than a certain level.

The background model will be estimated from a set of observations  $\chi_{ij} = \{x_{ij}(t-T), \dots, x_{ij}(t)\}$  where  $T$  (*bgs.nums*) is a time period used to adapt to changes. For each new sample, we update the training data set. According to this, new samples are added to the set and old ones are discarded, while the set size does not exceed a certain value.

We will model distribution of a particular pixel as a mixture of Gaussians following the technique proposed in<sup>3,4</sup>. Pixel values that do not fit the background distribution are considered foreground until there is a Gaussian that includes them with sufficient, evidence of supporting it. The  $M$  (*bgs.numg*) Gaussian mixture models can be expressed as:

$$p(x|\chi, BG+FG) = \sum_{m=1}^M \omega_m \mathcal{N}(\mu_m, \sigma_m^2 I), \quad (6)$$

where  $\mu_m$  and  $\sigma_m^2$  are the estimates of the mean and variance that describe the Gaussian component  $m$ . The covariance matrices are assumed to be diagonal and the identity matrix  $I$  has proper dimensions. Finally,  $\omega$  are positive mixing weights that add up to 1.

Given a new data sample  $x^t$  at the time  $t$ , the recursive update equations are:

$$\omega_m \leftarrow (1 - \alpha) \omega_m + \alpha o_m^t, \quad (7)$$

$$\mu_m \leftarrow (1 - p_m^t) \mu_m + p_m^t x^t, \quad (8)$$

$$\sigma^2_m \leftarrow (1 - p_m^t) \sigma^2_m + p_m^t \left( (x^t - \mu_m)^T (x^t - \mu_m) \right), \quad (9)$$

$$p_m^t = \frac{\alpha o_m^t}{\omega_m}, \quad (10)$$

where  $\alpha$  (*bgs.lstp*) describes a exponentially decaying envelope which is used to limit the influence of the old data and being approximately  $\alpha = 1/T$ . The ownership  $o_m^t$  is set to one for the “close” component with the largest  $\omega$ . We define that a sample is “close” to a component if the *Mahalanobis* distance from the component is, less than a particular value (*bgs.thre*). If there are no “close” components a new one is generated and if the maximum number of components is reached, we discard the component with smallest  $\omega$ .

Usually, the foreground objects will be represented by some additional clusters with small weights  $\omega$ . Therefore, we can approximate the background model by the first  $B$  largest clusters:

If the components are sorted to have descending weights, we have:

$$B = \arg \min_b \left( \sum_{m=1}^b \omega_m > (1 - c_f) \right), \quad (11)$$

where  $c_f$  is the maximum portion of the data that can belong to the foreground objects without influencing the background model.

## Animal detection

The algorithm to detect the animals in the tracking areas is defined as follows:

- Preprocessing: After the image has been obtained, the distortion map obtained in the calibration is used to create a distortion-free image of the arena. Then the background model removes the static parts of the image. The areas of the image outside the tracking areas defined by the user are also removed from the image.
- Segmentation: An intensity value (*det.thre*) is selected by the user to threshold the image into two binary sets. Finally, a closing mathematical morphological operation is executed. This operation removes small holes and imperfections detected bodies. The size of the structuring element (*det.elms*) and the iterations of the dilation (*det.dilt*) and erosion (*det.erot*).
- Filtering: The objects smaller or bigger than the minimum (*det.mins*) and maximum (*det.maxs*) size limits defined by the user are disregarded.

## Filtering

Additionally to the minimum (*det.mins*) and maximum (*det.maxs*) size limits defined by the user. The user can also define a number of additional operations and filters listed as follows:

- A closing morphological operation can be used setting the number of dilation operations (*det.dilt*), the number of erosion operations (*det.erot*), and the size of the structuring element used (*det.elss*).

- Objects can be filtered according to the radius of the corresponding minimum enclosing circle, setting a minimum (*det.minr*) and a maximum (*det.maxr*) size limit.
- Objects can be filtered according to the rate between the major and the minor radius of the minimum ellipse fitting the detected body, setting a minimum (*det.mish*) and a maximum (*det.mash*) size limits.
- Objects can be filtered according to the ratio between the area of the minimum ellipse fitting the body and the actual number of pixels detected in the body, setting a minimum threshold indicating the minimum fill rate of the object (*det.minf*).

## 2. Detection screen




Figure 11. Detection screen.

- 1:** Navigate through the frames of the video sequence.
- 2:** Displays previous frame.
- 3:** Displays next frame.
- 4:** Select previous sequence.
- 5:** Select next sequence.
- 6:** Select previous arena.
- 7:** Select next arena.
- 8:** Threshold, this parameter selects the minimum intensity level (1-255) for the detected objects.
- 9:** Selects minimum object size in pixels.
- 10:** Selects maximum object size in pixels.
- 11:** Allows to enable or disable the background subtraction algorithm (subtracted background will not be shown in the visualization window).
- 12:** Shows the original video frame.
- 13:** Shows the areas of the image excluded from the tracking area.
- 14:** Shows unfiltered detections.
- 15:** Shows filtered detections.
- 16:** Shows the number of detected objects with the current parameters.

**17:** Shows the range of sizes of detected objects, together with the mean and standard deviation of the sizes.

**18:** Displays current detection in the selected arena and video frame. Calibration distortion is applied to the displayed image. Detected objects are displayed in different colors over imposed to the image. The user should check that all animals are detected with the current parameters, and no objects of the background are detected.

## Tracking

### 1. Kalman tracking algorithm

Tracking is the problem of generating an inference about the motion of one or more objects from a sequence of images. The Kalman filter addresses the problem of estimating the state  $x \in R^n$  of a discrete-time controlled process that is governed by the linear equation, expressed as follows:

$$x(t+1) = Ax(t) + w(t), \quad (12)$$

where  $A$  is a  $n$  by  $n$  matrix called state transition matrix, which relates the state of the system at the previous time step to the state at the current step, and  $w$  represents the process noise, which is assumed normally distributed with mean 0.

For the state transition matrix, we consider the equations of two-dimensional motion assuming a constant acceleration between time steps:

$$x_{t+1} = x_t + v_{x,t} + \frac{1}{2} a_x, \quad (13)$$

$$y_{t+1} = y_t + v_{y,t} + \frac{1}{2} a_y, \quad (14)$$

$$v_{x,t+1} = v_{x,t} + a_x, \quad (15)$$

$$v_{y,t+1} = v_{y,t} + a_y, \quad (16)$$

where  $(x, y)$  is the animal position,  $(v_x, v_y)$  is the velocity and  $(a_x, a_y)$  is the acceleration, which is assumed constant in a time step. We also consider an observation model described by the following equation:

$$z(t) = Hx(t) + v(t), \quad (17)$$

where  $z \in R^m$  represents the measurement,  $H$  is a  $m$  by  $n$  matrix called observation matrix and  $v$  is the measurement error, which is assumed independent of  $w$  and normally distributed with mean 0.

So the model equations can be expressed as follows:

$$\begin{bmatrix} x_{t+1} \\ y_{t+1} \\ v_{x,t+1} \\ v_{y,t+1} \end{bmatrix} = \begin{bmatrix} 1 & 0 & t & 0 \\ 0 & 1 & 0 & t \\ 0 & 0 & 1 & 0 \\ 0 & 0 & 0 & 1 \end{bmatrix} \begin{bmatrix} x_t \\ y_t \\ v_{x,t} \\ v_{y,t} \end{bmatrix} + w(t), \quad (18)$$

$$\begin{bmatrix} z_x \\ z_y \end{bmatrix} = \begin{bmatrix} 1 & 0 & 0 & 0 \\ 0 & 1 & 0 & 0 \end{bmatrix} \begin{bmatrix} x_t \\ y_t \\ v_{x,t} \\ v_{y,t} \end{bmatrix} + v(t), \quad (19)$$

The Kalman filter works in a two-step recursive process. First, it estimates the new state, along with their uncertainties. Once the outcome of the next measurement (corrupted with noise) is observed, these estimates are updated. The algorithm can, run in real time using only the current input measurements and the previously calculated state. In the present work, the implementation of the Kalman filter was performed according to<sup>5</sup>, using an empirical estimate of the measurement error and the process noise covariances.

The Kalman filter essential problem is the assignment of detections to tracks. To this end, a cost is assigned to every possible pair of track–detection. The cost is understood as the probability of that detection to correspond to the current track position. It is calculated using the distance from the detected position to the predicted position of the animals. To this end, the minimum of the Euclidean distances is selected as cost metric according to the Hungarian optimization algorithm<sup>6</sup>.

To create a Kalman filter, we need to initialize the noise covariance matrices  $v$  and  $w$ , and the a posteriori error covariance matrix  $p$ . To this end we have used the following procedure.

$$w = \begin{bmatrix} dt^4 / 4 & 0 & dt^3 / 2 & 0 \\ 0 & dt^4 / 4 & 0 & dt^3 / 2 \\ dt^3 / 2 & 0 & dt^2 & 0 \\ 0 & dt^3 / 2 & 0 & dt^2 \end{bmatrix} a, \quad (20)$$

$$v = \begin{bmatrix} m & 0 \\ 0 & m \end{bmatrix}, \quad (21)$$

$$p = \begin{bmatrix} c & 0 & 0 & 0 \\ 0 & c & 0 & 0 \\ 0 & 0 & c & 0 \\ 0 & 0 & 0 & c \end{bmatrix}, \quad (22)$$

where  $dt$  (*kal.time*) represent a time increment magnitude,  $a$  (*kal.pron*) represent the estimated process noise,  $m$  (*kal.mean*) represent the estimated measurement noise, and  $c$  (*kal.errc*) is a value used to initialize the a posteriori error covariance matrix with a correct value. Modifying the  $a$  and  $m$  values will change the behavior of the filter, causing to use more aggressive predictions (not recommended) or to stick more to the measured positions. These values can be modified by user in the configuration text file.

## Feature Extraction

To be able to find the identity of individuals when multiple objects are tracked and there is risk of occlusion, we calculate a set of characteristic features for a detected body  $B$ . These features are the intensity histogram (HIST) and two Texture Center Maps (TCM) called the Intensity Center Map (ICM); and the Contrast Center Map (CCM), where the latter two are 2D distributions. They are defined as:

$$HIST_B(x) = \left| \{p : f_p = x\} \right|, \quad (23)$$

$$ICM_B(x, y) = \left| \{p : \|p - c\| = x, f_p + f_c = y\} \right|, \quad (24)$$

$$CCM_B(x, y) = \left| \{p : \|p - c\| = x, f_p + f_c = y\} \right|, \quad (25)$$

where  $p, c \in B$ , being  $p$  an arbitrary pixel and  $c$  the center of mass of the body  $B$ .  $f_p$  and  $f_c$  represent the color values of  $p$  and  $c$  respectively, and  $\|p - c\|$  is the Euclidean distance between  $p$  and  $c$ . Therefore, for a detected body in an image, we define the detection  $d$  as:

$$d(B) = \{c, |\{p\}|, HIST_B, ICM_B, CCM_B, T : p, c \in B\}, \quad (26)$$

where  $T$  is the time,  $|\{p\}|$  the size, and  $c$  is the position of the body.

The way these features are created and stored by the program can be altered. Therefore, the number of clusters of the histogram (*kal.hiss*), the maximum distance from the body center used by the TCM maps (*kal.tcmr*), the data type used to store the TCM maps in memory (*kal.tcnd*), and the maximum number of features stored for a track (*kal.hist*).

### Acceptance conditions

In order to reduce the possibility of assigning the wrong detection to a track, we have implemented an additional step to check if there is a significant change in the position or in the object size of the detections.

Assuming that the detection  $d_i$  has been assigned by the Hungarian algorithm to the track  $t_j$ , formed by the detections  $t_j = \{d_{j_0}, \dots, d_{j_n}\}$ , the following parameters are calculated:

**Frame distance (c1):** The distance from a detection  $d_i$  to a track  $f_i$ , is defined as the Euclidean distance from the predicted position of the track  $d_{j_{n+1}}$  to  $d_i$ .

$$c1 = \|t_j - d_i\| = \|d_{j_{n+1}} - d_i\|, \quad (27)$$

**Size change (c2):** The size of the detected body in  $d_i$  is compared with all detections in  $f_j$ , and the minimum relative change is selected

$$c2 = \max((|d_i| - |d_{j_n}|) / |d_{j_n}|, (|d_{j_n}| - |d_i|) / |d_i|), \quad (28)$$

To accept the assignment of  $d_i$  to  $t_j$  the following conditions must occur:

$$\begin{aligned} Accept &\leftarrow c2 < S_H \text{ and } c1 < d \times (T_i - T_{j_n}), \\ Accept &\leftarrow c2 < S_H \text{ and } c1 < d \end{aligned}, \quad (29)$$

where  $d$  (*kal.disf*) represent a frame distance condition,  $S_H$  (*kal.sich*) represent a size change condition.

## Collision

As explained in the main document, when two individuals overlap or cross. We mark the conflicted tracks as inactive tracks and generate new ones. To detect a collision we define the operator  $cl(t_j) = d_i$  to refer to the closest detection  $d_i$  to the track  $t_j$ .

A track is classified as conflicted if one of the following conditions applies:

$$cl(t_i) = cl(t_k) = d_i \text{ and } abs(\|t_k - d_i\| - \|t_k - d_i\|) < (d \times a_d), \quad (30)$$

$$abs(\|t_k - d_i\| - \|t_k - d_i\|) < a, \quad (31)$$

where  $a_d$  (*kal.advr*) is defined from 0 to 1, and  $a$  (*kal.advm*) represent a minimum advantage value in pixels.

## Delete conditions

To avoid that occasional missdetections interfere with the final results, tracks that accomplish one of the following conditions are deleted.

- Tracks which are non-active and are smaller than the minimum age (*kal.dmax*). If the size of the track is smaller than (*kal.mins*) the track is also marked as short.
- Tracks not assigned to any detection for a certain amount of frames (*kal.dage*) are marked as inactive. If the size of the track is smaller than (*kal.mins*) the track is also marked as short.

## 2. Fragment Identification

This is an optional post-processing step that calculates what trajectory fragments that belong to each individual, preserving the identity of the animals after an occlusion. To assign the correct identities, we compute and study the similarities of the trajectory features. First, an identity matrix is constructed, containing a similarity value for each pair of tracks  $Sim(t_{row}, t_j)$ . This value is constructed comparing the stored features of those tracks, and represent the likelihood of those tracks to correspond to the same animal.

Since to estimate the similarity of two tracks is a very computationally expensive procedure, only the samples similar in size (*kal.cmse*), and with a minimum histogram correlation value (*kal.cmhc*) are used to construct the similarity values  $Sim(t_{row}, t_j)$

To assign the different tracks to each other, the first step is to select groups of long tracks, coexisting at the same time and belonging to all individuals. We know this is an assignment problem with an optimal solution, so we use a variant of the classic Hungarian optimization algorithm (Kuhn 1955). The minimum similarity value (*kal.idgb*) accepted by the algorithm can be modified by the user.

To assign the remaining tracks, we iteratively select the best correlation value in the *matrix*, first with the long and then with the short tracks. The minimum similarity values (*kal.idlb*, *kal.idsrb*) and the minimum average similarity values (*kal.idla*, *kal.idsa*) accepted can also be modified by the user.

With every assignment, we update the matrix propagating the knowledge obtained iteration and reducing the uncertainty for the remaining tracks.

Since the fragment identification algorithm works as a post processing technique, and requires a big amount of memory to store the features of active and inactive tracks, the number of tracks kept in memory can be defined by the user (*kal.idff*). When the number of tracks reaches this number, the algorithm will try to identify the tracks to release memory, before continuing the tracking

### 3. Tracking screen




Figure 12. Tracking Screen.

- 1: Select the number of animals in each arena. (All arenas should have the same number of animals, or be empty).
- 2: Selects fragment identification algorithm. Available options are: “**Not Id**” (no fragment identification algorithm will be used) “**Hist. sel. by Shape (beta)**” (only histogram information will be used in the fragment identification algorithm. Reduces processing time and memory use significantly but this algorithm is still under development), and “**2TCM sel. by Hist.**” (uses the fragment identification algorithm proposed in the paper).
- 3: Selects the maximum pixels per frame, an animal is allowed to move. It is recommended to overestimate this parameter.
- 4: Enables the user to change advanced options (5-16).
- 5: Use Avg. Use a mean metric instead of a max similarity value (Recommended).
- 6: Use High order Correlation. Takes in account the similarity with all fragments, without reducing the speed. (Recommended).
- 7: Allow all groups to be matched. Increases accuracy, but may lead to have a final number of tracks different to the number of individuals.
- 8: History, maximum number of samples used for every track in the algorithm. This parameter also affects significantly the memory and time used by the fragment identification algorithm.
- 9: Comparisons, limits the number of comparisons between two samples of the same track in order to increase speed.
- 10: Size (number of color clusters) used by the histograms.
- 11: Max distance from the body center used by the TCM maps. Setting a higher the value, will allow to use a greater part of the animal bodies in the algorithm, increasing accuracy, memory use and required time.
- 12: Minimum track size to be used by the algorithm, smaller tracks will not be saved.
- 13: Minimum track size to classify it as long.

## Results

### 1. Statistics screen

#### Individual statistics




Figure 13. Statistics screen, individuals.

- 1: Change between information of one single individual and for the entire population.
- 2: Select previous sequence.
- 3: Select next sequence.
- 4: Select previous arena.
- 5: Select next arena.
- 6: Displays the tracking statistics for the selected arena and sequence.
- 7: Enables the user to change advanced options (8-11).
- 8: Allows the user to change the speed threshold to estimate the mobility rate of the animal in the stats.
- 9: Allow the user to change the distance threshold (in mm) to calculate the frozen events. A frozen event is detected when the animal has moved less than this value in time interval.
- 10: Allow the user to change the time threshold (in seconds) to calculate the frozen events. A frozen event is detected when the animal has moved less than a distance threshold during a lapse of time.
- 11: Allow the user to change the time threshold (in seconds) to calculate a transition. A transition is detected when the time between two consecutive detections exceed this value, signaling a period of time when the animal is not visible.
- 12: Saves all results in the output folder.
- 13: Shows a free resizable window showing the tracking for the current arena, the window can be closed pressing the 'esc' key, or clicking in the window upper right corner.
- 14: If enabled the generate output button will also save a video of the tracking in the output folder.
- 15: Select the starting point (for all sequences) to show the tracking.
- 16: Select the ending point (for all sequences) to show the tracking.

## Population statistics




Figure 14. Statistics screen, population.

**1:** Displays the tracking statistics for the entire population.

## 2. Graphical outputs screen

### Individual graphical outputs




Figure 15. Graphical outputs, individuals.

**1:** Selects one of the different graphic available results. The graphical outputs will be shown superimposed to the arena picture, and in real scale. The heat maps are color coded according to the normalized representation of the frequency of use, in the linear scale shown in Figure 19.

Available options are:

- “**Plain Output**”. This option will show the first frame of the sequence and arena without any graphical output.
- “**North Edge**”. This option will show a heat map showing the use of different zones in the arena according to the distance of to the North (N) wall. See section: “Distance to Edges (*Dist\_Edges.txt, Dist\_Edge\_[edge].jpeg*)”.
- “**West Edge**”. This option will show a heat map showing the use of different zones in the arena according to the distance of to the North (W) wall. See section: “Distance to Edges (*Dist\_Edges.txt, Dist\_Edge\_[edge].jpeg*)”.
- “**South Edge**”. This option will show a heat map showing the use of different zones in the arena according to the distance to the North (S) wall. See section: “Distance to Edges (*Dist\_Edges.txt, Dist\_Edge\_[edge].jpeg*)”.
- “**East Edge**”. This option will show a heat map showing the use of different zones in the arena according to the distance to the North (E) wall. See section: “Distance to Edges (*Dist\_Edges.txt, Dist\_Edge\_[edge].jpeg*)”.
- “**All Edges**”. This option will show a heat map showing the use of different zones in the arena according to the distance to any wall. See section: “Distance to Edges (*Dist\_Edges.txt, Dist\_Edge\_[edge].jpeg*)”.
- “**Exploration**”. This option shows a heat map showing the use of different zones in the arena using a regular grid with a grid of a size selected by the user. See section: “Exploration (*Exploration.txt, Exploration.jpeg*)”.
- “**Center Point Dist.**”. This option will show a heat map showing the use of different zones in the arena according to the distance to the center of the arena. See section: Distance to Center Position (*Dist\_CenterPos.txt, Dist\_Center\_Pos.jpeg*).
- “**Mean Point Dist.**”. This option will show a heat map showing the use of different zones in the arena according to the distance to the mean point of the detected animals in the arena. See section: “Distance to Mean Position (*Dist\_MeanPos.txt, Dist\_Mean\_Pos.jpeg*)”.
- “**Trajectory**”. This option shows a representation of the trajectory of the animals. See section: “Tracking in Real Space coordinates (*Tracking\_RealSpace.txt, Trajectory.jpeg*)”.

**2:** Change zone size in mm. This will affect the zone statistics.

**3:** Number of zones. This will affect the edge zone statistics.

**4:** Normalize arenas. If enabled, the zones will be computed according to the extremes of the detected positions instead of using the full image.

**5:** Visualization panel, displays the current graphic selected in the first image of the current video, for the current arena.

## Population graphical outputs




Figure 16. Graphical outputs, population.

**1:** If the arenas are not symmetrical, and do not have the same orientation. This button allows to change the orientation of the arenas, in the projected virtual arena for the entire population. Available options are: “**Same Orientation**”, “**Vertical Mirror**” (the orientation of the arenas in the upper and lower half of the image is inverted in its vertical axis), “**Horizontal Mirror**” (the orientation of the arenas in the left and right half of the image is inverted in its horizontal axis), “**Vertical & Horizontal Mirror**” (combines the previous two options).

**2:** Visualization panel, displays the selected graphic in a real scale virtual arena, projecting together the results of the entire population.

### 3. Output Files

The application automatically generates a folder structure to save the results of the tracking. A set of files are generated to save the project configuration data (Table 3). The main results of the project, containing the majority of the statistical data generated for each individual and the population will be formatted in a datasheet (Figure 4). Additionally, a set of files is created with the results for the entire population (Table 5) and another set of files is created with the results for each arena and sequence (Table 6). The creation of most of these files can be enabled or disabled in the *Configuration.txt* file. When analyzing a project, a folder is generated for every video sequence and the corresponding file names are named ending in a number representing the arena number.

Table 3. Structure of the output project files, these files are detailed in Table 2.

Output Folder	Output Files
C:/ProjectFolder/	pname.tox
C:/ProjectFolder/pname/	pname_Input.txt pname_Configuration.txt pname_Arena.txt pname_ArenaNames.txt pname_Calibrator.txt pname_Output.txt

Table 4. Datasheet continuing main results of the experiment.

File Name	Video	Seq	Output Folder	Output Excel File
Seq1_0.avi	1	1	C:/ProjectFolder/	pname.xls
Seq1_1.avi	2			
Seq2_0.avi	1			

Table 5. Structure of the output files containing the results for the entire population.

File Name	Video	Seq	Output Folder	Output TEXT Files	Output JPG Files
Seq1_0.avi	1	1	C:/ProjectFolder/pname	Tracking.txt Tracking_RealSpace.txt Instant_Accel.txt Instant_Speed.txt Dist_MeanPos.txt Dist_CenterPos.txt Dist_Edges.txt Exploration.txt Transitions.txt FrozenEvents.txt Stats.txt	Dist_CenterPos.jpeg Dist_Edge_All.jpeg Dist_Edge_E.jpeg Dist_Edge_N.jpeg Dist_Edge_S.jpeg Dist_Edge_W.jpeg Dist_MeanPos.jpeg Exploration.jpeg Trajectory.jpeg
Seq1_1.avi	2				
Seq2_0.avi	1				

Table 6. Structure of the output files containing the results for each arena.

File Name	Video	Seq	Output Folder	Arena	Output TEXT Files	Output JPEG Files
Seq1_0.avi	1	1	C:/ProjectFolder/pname /Seq1/	1	Tracking_0.txt Tracking_RealSpace_1.txt Instant_Accel_1.txt Instant_Speed_1.txt Dist_MeanPos_1.txt Dist_CenterPos_1.txt Dist_Edges_1.txt Exploration_1.txt Transitions_1.txt FrozenEvents_1.txt	Dist_CenterPos.jpeg Dist_Edge_All.jpeg Dist_Edge_E.jpeg Dist_Edge_N.jpeg Dist_Edge_S.jpeg Dist_Edge_W.jpeg
Seq1_1.avi	2				Tracking_1.txt Tracking_RealSpace_2.txt Instant_Accel_2.txt Instant_Speed_2.txt Dist_MeanPos_2.txt Dist_CenterPos_2.txt Dist_Edges_2.txt Exploration_2.txt Transitions_2.txt FrozenEvents_2.txt	Dist_MeanPos.jpeg Exploration.jpeg Trajectory.jpeg

### Tracking (*Traking.txt*)

This file is used by the program. And contains the animal detected positions. Coordinates are in pixel and relative to the individual after removing distortion. So they cannot be back-projected directly to the original images. Time coordinates are in frames. An example and the meaning of the columns is the following. The arenas are numbered from 0 to n-1, in accordance to the c++ vector indexing.

Table 7. *Traking.txt* files.

Frame Number	Arena Number	Track number	X-Position	Y-Position
360	0	1	224.634	32.2103
361	0	1	225.131	32.2103
362	0	1	226.737	32.2103
363	0	1	227.434	32.2103
364	0	1	226.624	32.2103

Label	Meaning
0	Predicted Position
1	Confirmed position
2	Occluded position
3	Mirror position

### Tracking in Real Space coordinates (*Tracking\_RealSpace.txt*, *Trajectory.jpeg*)

This file is the direct translation of the previous file to a real space. Coordinates are expressed in millimeters (typically, though it depends on the calibration unit used) and time in seconds.

In the files named *Tracking\_RealSpace\_[arena\_number].txt* space coordinates will not be referred to the arena, but to a reference point depending on the pose, so different arenas will have different spatial coordinates, corresponding to their relative position in the real space.

In the file related to the population, *Tracking\_RealSpace.txt* the space coordinates will be projected to a common virtual arena, according to a selected arena distribution. If the normalized option is selected (*ana.norm*) so the min and x and y coordinates from each arena will correspond to the same x and y coordinates from the virtual arena. The arenas are numbered from 0 to n-1Here, and in consequent results, the arenas are numbered from 1 to n, in accordance with the values shown in the interface.

Table 8. *Tracking\_RealSpace.txt* files.

Time (sec)	Arena	Track	Pos. X (mm)	Pos. Y (mm)	Label
15.0414	1	1	65.2423	31.7268	1
15.0414	1	1	64.8647	31.7615	1
15.0831	1	1	64.1963	31.7955	1
15.1248	1	1	63.6246	31.7977	1
15.1664	1	1	63.155	31.7606	1
15.2081	1	1	62.7215	31.6936	1
15.2498	1	1	62.3192	31.6089	1
15.2914	1	1	61.9783	31.5102	1
15.3331	1	1	61.7018	31.4046	1
15.3748	1	1	61.4741	31.2961	1

The graphic outputs, *Trajectory.jpeg* correspond to colored representation of the trajectory of the animal, superimposed to the arena picture, and in real scale. The trajectory of different tracks in the same arena is colored with a different color.




Figure 17. *Trajectory.jpeg*. Trajectory projection of 5 fish in a single arena experimental setup.




Figure 18. *Trajectory.jpeg*. Trajectory projection of 3 fish in a four arena experimental setup.

### Instantaneous Speed (*Instant\_Speed.txt*)

This file is the first derivate of the *Tracking\_RealSpace.txt* files. Where the instantaneous speed  $s_i$  corresponding to the time  $t_i$  is calculated according to the formula:

$$s_i^{arena,track} = \frac{\sqrt{(x_{i+c} - x_{i-c})^2 + (y_{i+c} - y_{i-c})^2}}{t_{i+c} - t_{i-c}}, \quad (32)$$

where  $c$  (*ana.spsa*) is a sampling distance defined by the user which can be tuned by the user.

Table 9. *Instant\_Speed.txt* files.

Time (sec)	Arena	Track	Current Speed (mm/sec)
0.08	1	1	73.2125
0.12	1	1	59.5829
0.16	1	1	52.8199
0.2	1	1	43.0816
0.24	1	1	32.5418
0.28	1	1	23.6301
0.32	1	1	17.503
0.36	1	1	13.6693
0.4	1	1	10.7872
0.44	1	1	13.3477

### Instantaneous Acceleration (*Instant\_Accel.txt*)

This file is the second derivate of the *Tracking\_RealSpace.txt* files. Where the instantaneous acceleration  $a_i$  corresponding to the time  $t_i$  is calculated according to the formula (accelerations and decelerations will be represented as positive values):

$$a_i^{arena,track} = \frac{\sqrt{(s_{i+c} - s_{i-c})^2}}{t_{i+c} - t_{i-c}}, \quad (33)$$

where  $c$  (*ana.spsa*) is a sampling distance defined by the user which can be tuned by the user.

Table 10. *Instant\_Accel.txt* files.

Time (sec)	Arena	Track	Current Accel. (mm/s^2)
0.16	1	1	254.192
0.2	1	1	224.705
0.24	1	1	220.731
0.28	1	1	183.827
0.32	1	1	135.966
0.36	1	1	64.265
0.4	1	1	39.7648
0.16	1	1	254.192
0.2	1	1	224.705
0.24	1	1	220.731

### Distance to Edges (*Dist\_Edges.txt*, *Dist\_Edge\_[edge].jpeg*)

We define the walls of the arena by the lines connecting its corner points. To this end we will use a cardinal nomenclature. Therefore the North (N) wall will be defined as the line connecting the North West (NW) and North East (NE) corners of the image and it will correspond to the up side of the arena image.

Therefore, the distance for a detected position  $(x_i, y_i)$  corresponding to the time  $t_i$ , the distance corresponding to the different walls of the arena will be defined as follows:

$$dist_{N,i} = \frac{|(y_{NE} - y_{NW})x_i - (x_{NE} - x_{NW})y_i + x_{NE}y_{NW} - y_{NE}x_{NW}|}{\sqrt{(x_{NE} - x_{NW})^2 + (y_{NE} - y_{NW})^2}}, \quad (34)$$

$$dist_{W,i} = \frac{|(y_{SW} - y_{NW})x_i - (x_{SW} - x_{NW})y_i + x_{SW}y_{NW} - y_{SW}x_{NW}|}{\sqrt{(x_{SW} - x_{NW})^2 + (y_{SW} - y_{NW})^2}}, \quad (35)$$

$$dist_{S,i} = \frac{|(y_{SE} - y_{SW})x_i - (x_{SE} - x_{SW})y_i + x_{SE}y_{SW} - y_{SE}x_{SW}|}{\sqrt{(x_{SE} - x_{SW})^2 + (y_{SE} - y_{SW})^2}}, \quad (36)$$

$$dist_{E,i} = \frac{|(y_{NE} - y_{SE})x_i - (x_{NE} - x_{SE})y_i + x_{NE}y_{SE} - y_{NE}x_{SE}|}{\sqrt{(x_{NE} - x_{SE})^2 + (y_{NE} - y_{SE})^2}}, \quad (37)$$

We will then assign the detection  $i$  to an area  $k$  for each wall (and for any wall), according to the distance of the point to the wall as defined as follows:

$$N : i \in k \quad if \quad k(dst) < dist_{N,i} \leq k(dst+1), \quad (38)$$

$$W : i \in k \quad if \quad k(dst) < dist_{W,i} \leq k(dst+1), \quad (39)$$

$$E : i \in k \quad if \quad k(dst) < dist_{E,i} \leq k(dst+1), \quad (40)$$

$$\text{All} : i \in k \quad if \quad k(dst) < \min(\{dist_{N,i}, dist_{W,i}, dist_{S,i}, dist_{E,i}\}) \leq k(dst+1), \quad (41)$$

where  $dst$  (*ana.zsiz*) is an arbitrary distance given by the user to define the number of areas. Therefore, for every arena, we count all the assignments for each area and present the results in three different metrics, the total number of detections (frame count, Table 11), the total number of detections divided by the frame rate (representing the amount of seconds in each area, Table 12), and the frequency in each area (the number of detections in that area divided by the total number of detections, Table 13).

Table 11. *Dist\_Edges.txt* files, total number of detections.

	Frame Count									
	20mm	40mm	60mm	80mm	100mm	120mm	140mm	160mm	180mm	200mm
Edge N	205850	159661	112987	89949	77866	72482	71472	67567	67440	70892
Edge W	182286	150407	115214	86531	75916	71105	67611	66397	68303	70350
Edge S	191830	146005	100341	81174	71692	67630	67383	70903	72280	77001
Edge E	212245	152288	96947	81260	73388	68878	66145	67615	69283	73340
Edge ALL	747518	412348	182703	91367	51580	28616	11324	118	0	0

Table 12. *Dist\_Edges.txt* files, total number of detections divided by the frame rate.

	Time Count (sec)									
	20mm	40mm	60mm	80mm	100mm	120mm	140mm	160mm	180mm	200mm
Edge N	8234	6386.44	4519.48	3597.96	3114.64	2899.28	2858.88	2702.68	2697.6	2835.68
Edge W	7291.44	6016.28	4608.56	3461.24	3036.64	2844.2	2704.44	2655.88	2732.12	2814
Edge S	7673.2	5840.2	4013.64	3246.96	2867.68	2705.2	2695.32	2836.12	2891.2	3080.04
Edge E	8489.8	6091.52	3877.88	3250.4	2935.52	2755.12	2645.8	2704.6	2771.32	2933.6
Edge ALL	29900.7	16493.9	7308.12	3654.68	2063.2	1144.64	452.96	4.72	0	0

Table 13. *Dist\_Edges.txt* files, frequency in each area (the number of detections in that area divided by the total number of detections).

	Frequency (Zone Detections/Total Detections)									
	20mm	40mm	60mm	80mm	100mm	120mm	140mm	160mm	180mm	200mm
Edge N	0.13	0.10	0.07	0.06	0.05	0.05	0.05	0.04	0.04	0.05
Edge W	0.12	0.10	0.08	0.06	0.05	0.05	0.04	0.04	0.04	0.05
Edge S	0.13	0.10	0.07	0.05	0.05	0.04	0.04	0.05	0.05	0.05
Edge E	0.14	0.10	0.06	0.05	0.05	0.05	0.04	0.04	0.05	0.05
Edge ALL	0.49	0.27	0.12	0.06	0.03	0.02	0.01	0.00	0.00	0.00

Note: In practice, for memory reason the matrix size will be predefined, and the output can be bigger than the actual arena being filled with zeros.

The graphic outputs, *Dist\_Edge\_[edge].jpeg* correspond to colored representation of the frequency of the use of each zone. To this end the matrix is normalized so the maximum frequency will be represented as 1 and the minimum none zero value will be coded as a close to zero constant. The color scale is linear and from low to high frequency, the color codes used will go from blue to green to yellow to red as shown in Figure 18. In the graphic outputs, the zones will be superimposed to the images of the arenas, and the zone size is represented in a true scale.




Figure 19. Color linear scale for representing the normalized of the frequency of use of a zone in the heat maps.




Figure 20. *Dist\_Edge\_[edge].jpeg*. Representation of the zone use as heat maps according to de distance to the edges.

### Exploration (*Exploration.txt, Exploration.jpeg*)

The arena is divided in regular square non overlapping zones of a size selected by user. And the use of each zone is computed for each arena in the same way as the edge zones. Therefore a matrix of zones will be created, so the zone k corresponding to the coordinates of the zone matrix will be assigned to a detected position  $(x_i, y_i)$  corresponding to the time  $t_i$  according to the following formula:

$$i \in k(k_x, k_y) \left\{ \begin{array}{l} \text{floor} \left( \frac{\left| (y_{NE} - y_{NW})x_i - (x_{NE} - x_{NW})y_i + x_{NE}y_{NW} - y_{NE}x_{NW} \right|}{\sqrt{(x_{NE} - x_{NW})^2 + (y_{NE} - y_{NW})^2}} \right) \\ \text{floor} \left( \frac{\left| (y_{SW} - y_{NW})x_i - (x_{SW} - x_{NW})y_i + x_{SW}y_{NW} - y_{SW}x_{NW} \right|}{\sqrt{(x_{SW} - x_{NW})^2 + (y_{SW} - y_{NW})^2}} \right) \end{array} \right\}, \quad (42)$$

where  $dst$  (*ana.zsiz*) represents the size of each square of the arena, selected by the user.

For every arena, we count all the assignments for each area and present the results in three different metrics, the total number of detections (frame count, Table 14), the total number of detections divided by the frame rate (representing the amount of seconds in each area, Table 15), and the frequency in each area (the number of detections in that area divided by the total number of detections, Table 16).

Table 14. *Exploration.txt* files, total number of detections.

	Frame Count									
	20mm	40mm	60mm	80mm	100mm	120mm	140mm	160mm	180mm	200mm
20mm	10865	16504	15511	13738	13299	13319	12400	12151	13146	14016
40mm	17455	14569	11358	10022	8670	7977	7707	7670	7777	8864
60mm	16034	11226	12950	5912	5334	4370	4135	4164	4206	4057
80mm	15243	9955	6676	4223	3629	3317	2928	2627	2786	2918
100mm	13531	8535	5285	3121	2752	2751	2508	2292	2479	2507
120mm	12313	8443	5249	2789	2708	2641	2274	2261	2081	2176
140mm	11653	8242	4251	3125	2214	2292	2496	2268	1996	2177
160mm	10503	8381	4527	3424	2234	2196	2301	2263	1884	1912
180mm	10561	7757	4627	3391	2577	2242	2234	2158	2084	2039
200mm	11596	7637	5031	3309	2963	2838	2228	2140	2383	2299

Table 15. *Exploration.txt* files, total number of detections divided by the frame rate.

	Time Count (sec)									
	20mm	40mm	60mm	80mm	100mm	120mm	140mm	160mm	180mm	200mm
20mm	434.6	660.2	620.4	549.5	532.0	532.8	496.0	486.0	525.8	560.6
40mm	698.2	582.8	454.3	400.9	346.8	319.1	308.3	306.8	311.1	354.6
60mm	641.4	449.0	518.0	236.5	213.4	174.8	165.4	166.6	168.2	162.3
80mm	609.7	398.2	267.0	168.9	145.2	132.7	117.1	105.1	111.4	116.7
100mm	541.2	341.4	211.4	124.8	110.1	110.0	100.3	91.7	99.2	100.3
120mm	492.5	337.7	210.0	111.6	108.3	105.6	91.0	90.4	83.2	87.0
140mm	466.1	329.7	170.0	125.0	88.6	91.7	99.8	90.7	79.8	87.1
160mm	420.1	335.2	181.1	137.0	89.4	87.8	92.0	90.5	75.4	76.5
180mm	422.4	310.3	185.1	135.6	103.1	89.7	89.4	86.3	83.4	81.6
200mm	463.8	305.5	201.2	132.4	118.5	113.5	89.1	85.6	95.3	92.0

Table 16. *Exploration.txt* files, frequency in each area (the number of detections in that area divided by the total number of detections).

	Frequency (Zone Detections/Total Detections)									
	20mm	40mm	60mm	80mm	100mm	120mm	140mm	160mm	180mm	200mm
20mm	0.0071	0.0108	0.0102	0.0090	0.0087	0.0087	0.0081	0.0080	0.0086	0.0092
40mm	0.0114	0.0095	0.0074	0.0066	0.0057	0.0052	0.0051	0.0050	0.0051	0.0058
60mm	0.0105	0.0074	0.0085	0.0039	0.0035	0.0029	0.0027	0.0027	0.0028	0.0027
80mm	0.0100	0.0065	0.0044	0.0028	0.0024	0.0022	0.0019	0.0017	0.0018	0.0019
100mm	0.0089	0.0056	0.0035	0.0020	0.0018	0.0018	0.0016	0.0015	0.0016	0.0016
120mm	0.0081	0.0055	0.0034	0.0018	0.0018	0.0017	0.0015	0.0015	0.0014	0.0014
140mm	0.0076	0.0054	0.0028	0.0020	0.0015	0.0015	0.0016	0.0015	0.0013	0.0014
160mm	0.0069	0.0055	0.0030	0.0022	0.0015	0.0014	0.0015	0.0015	0.0012	0.0013
180mm	0.0069	0.0051	0.0030	0.0022	0.0017	0.0015	0.0015	0.0014	0.0014	0.0013
200mm	0.0076	0.0050	0.0033	0.0022	0.0019	0.0019	0.0015	0.0014	0.0016	0.0015

The graphic output, *Exploration.jpeg* represent the frequency of use of each area in the same way as explained before.




Figure 21. *Exploration.jpeg*. Representation of the zone use (exploration) as a heat map with a grid of a size selected by the user, for a single arena experimental setup with 5 fish.




Figure 22. *Exploration.jpeg*. Representation of the zone use (exploration) as a heat map with a grid of a size selected by the user for a four arena experimental setup with three fish.

### Distance to Mean Position (*Dist\_MeanPos.txt*, *Dist\_Mean\_Pos.jpeg*)

This output is similar to *Dist\_Edges.txt*, representations, with the difference that the zones will be created according to the distance of a detected position  $(x_i, y_i)$  corresponding to the time  $t_i$ , to the mean of the detected positions.  $(x_M, y_M)$ . This is defined as follows.

$$dist_{M,i} = \sqrt{(x_M - x_i)^2 + (y_M - y_i)^2}, \quad (43)$$

We will assign the detection  $i$  to an area  $k$ , according to the distance of the point to the mean position defined as follows:

$$M : i \in k \quad \text{if} \quad k(dst) < dist_{M,i} \leq k(dst + 1), \quad (44)$$

where  $dst$  (*ana.zsizq*) is an arbitrary distance given by the user to define the number of areas.

For every arena, we count all the assignments for each area and present the results in three different metrics, the total number of detections (frame count, Table 17), the total number of detections divided by the frame rate (representing the amount of seconds in each area, Table 18), and the frequency in each area (the number of detections in that area divided by the total number of detections, Table 19).

Table 17. *Dist\_MeanPos.txt* files, total number of detections.

	20mm	40mm	60mm	80mm	100mm	120mm	140mm	160mm	180mm	200mm
Frame Count	7264	20741	35441	59533	101590	200553	445979	371775	240955	41882

Table 18. *Dist\_MeanPos.txt* files, total number of detections divided by the frame rate.

	20mm	40mm	60mm	80mm	100mm	120mm	140mm	160mm	180mm	200mm
Time Count (sec)	290.56	829.64	1417.64	2381.32	4063.6	8022.12	17839.2	14871	9638.2	1675.28

Table 19. *Dist\_MeanPos.txt* files, frequency in each area (the number of detections in that area divided by the total number of detections).

	20mm	40mm	60mm	80mm	100mm	120mm	140mm	160mm	180mm	200mm
Frequency (Zone Det./Total Det.)	0.0048	0.0136	0.0232	0.0390	0.0666	0.1314	0.2923	0.2437	0.1579	0.0275

The graphic output, *Dist\_Mean\_Pos.jpeg* represent the frequency of use of each area in the same way as explained before.




Figure 23. *Dist\_Mean\_Pos.jpeg*. Representation of the zone use as heat maps according to de distance to mean position of the arena, for a single arena experimental setup with 5 fish.




Figure 24. *Dist\_Mean\_Pos.jpeg*. Representation of the zone use as heat maps according to de distance to mean position of the arena, for a four arena experimental setup with three fish.

### Distance to Center Position (*Dist\_CenterPos.txt*, *Dist\_Center\_Pos.jpeg*)

This output is similar to the *Dist\_Edges.txt*, representations, with the difference that the zones will be created according to the distance of a detected position  $(x_i, y_i)$  corresponding to the time  $t_i$ , to the center of the arena.  $(x_C, y_C)$ . This is defined as follows.

$$dist_{C,i} = \sqrt{(x_C - x_i)^2 + (y_C - y_i)^2}, \quad (45)$$

We will assign the detection  $i$  to an area  $k$ , according to the distance of the point to center position defined as follows:

$$C : i \in k \quad \text{if} \quad k(dst) < dist_{C,i} \leq k(dst + 1). \quad (46)$$

where  $dst$  (*ana.zsizq*) is an arbitrary distance given by the user to define the number of areas.

For every arena, we count all the assignments for each area and present the results in three different metrics, the total number of detections (frame count, Table 10), the total number of detections divided by the frame rate (representing the amount of seconds in each area, Table 21), and the frequency in each area (the number of detections in that area divided by the total number of detections, Table 21).

Table 20. *Dist\_CenterPos.txt* files, total number of detections.

	20mm	40mm	60mm	80mm	100mm	120mm	140mm	160mm	180mm	200mm
Frame Count	7226	20643	35576	59870	100260	199441	445106	373978	244069	39543

Table 21. *Dist\_CenterPos.txt* files, total number of detections divided by the frame rate.

	20mm	40mm	60mm	80mm	100mm	120mm	140mm	160mm	180mm	200mm
Time Count (sec)	289.04	825.72	1423.04	2394.8	4010.4	7977.64	17804.2	14959.1	9762.76	1581.72

Table 22. *Dist\_CenterPos.txt* files, frequency in each area (the number of detections in that area divided by the total number of detections).

	20mm	40mm	60mm	80mm	100mm	120mm	140mm	160mm	180mm	200mm
Frequency (Zone Det./Total Det.)	0.0047	0.0135	0.0233	0.0392	0.0657	0.1307	0.2917	0.2451	0.1600	0.0259

The graphic output, *Dist\_Center\_Pos.jpeg* represent the frequency of use of each area in the same way as explained before.




Figure 25. *Dist\_Center\_Pos.jpeg*. Representation of the zone use as heat maps according to de distance to central position of the arena, for a single arena experimental setup with 5 fish.




Figure 26. *Dist\_Center\_Pos.jpeg*. Representation of the zone use as heat maps according to de distance to central position of the arena, for a four arena experimental setup with three fish.

### Transitions (*Transitions.txt*)

Transitions, count the moments when tracked objects appears in the image (transitions labeled as 1), or disappear from it (transitions labeled as 0).

The time, track, position and other parameters from the first/last detection is saved in a file.

A transition is detected when the time between two consecutive detections exceed a value (*ana.ttim*) selected by the user (for the first and last detection the time is computed from the video start or the video end time).

Table 23. *Transitions.txt* files.

Time (sec)	Video Seq.	Arena	Track	Pos. X (mm)	Pos. Y (mm)	Label
1.92	0	1	1	30.1997	139.79	0
16.16	0	1	2	40.3914	142.042	1
25.56	0	1	2	264.587	142.659	0
36.28	0	1	3	67.851	142.557	1
43.92	0	1	3	260.846	141.895	0
59.8	0	1	4	37.3171	142.579	1
63	0	1	4	26.5074	141.028	0
75.44	0	1	5	266.598	142.652	1
78.32	0	1	5	268.247	141.085	0
86.12	0	1	6	42.6465	142.008	1

### Frozen Events (*FrozenEvents.txt*)

A frozen event represent the moments when the animal stands still in a predetermined position of the arena for a period of time.

The time, track number, position, length and other parameters from event are saved in a file.

A frozen event is detected when the animal has moved less than a value (*ana.fnmmt*) in a certain amount of time (*ana.ftim*) selected by the user

Table 24. *FrozenEvents.txt* files.

Time (sec)	Video Seq.	Arena	Track	Avg. Pos. X (mm)	Avg. Pos. Y (mm)	Time Length (sec)
12.16	0	1	1	55.0868	33.9619	3.52
98.68	0	1	1	105.536	41.5467	3.6
37.12	6	1	2	188.977	74.738	7.32
51.92	6	1	2	195.399	76.9609	8.88
65.96	6	1	2	199.721	75.0917	4.6
79.32	6	1	2	208.338	75.3267	4
5.8	11	1	3	235.564	28.7672	8.24
21.16	11	1	3	223.176	27.4975	5
29.64	11	1	3	214.962	27.6328	5.8
39.08	11	1	3	202.831	27.8868	13.2

### Stats (*Stats.txt*)

The *Stats.txt* files represent a summary of all the relevant parameters of the arena (Table 25), and summarize the tracking results with global statistics (Table 26). For the population the main and standard deviation of parameters are calculated.

Table 25. *Stats.txt* files, sequence and arena information and video data.

Parameter Name	Example of Value	Description
Video Resolution	[2048 x 2048]	Size in pixels of a video frame
Video FrameRate	25	Number of frames per second in the original video
Analysed Video Frames	22500	Number of processed video frames (assumed to be equal for every sequence)
Analysed Video Time	900	(Analysed Video Frames)/ (Video Framerate)
Arena	1 Arena1 [0, 0] [990.332, 0], [0, 991.625] [990.332, 991.625] [0, 0, 1] [287.538, 0, 1] [0, 282.707, 1] [287.538, 282.707, 1]	Arena number and arena name. Corners delimiting the arena in 2D space and in the 3D World, in this order: [NW] [NE] [SW] [SE]
Arena Size	[1506 x 1078] [1506 x 1078]	The dimensions of the arena in pixels and in real scale units.
Arena Center	[141.354, 143.769]	Position in mm of the center of the arena
Mean Position	[144.571, 4.32984]	Position in mm of the mean of all detections

Table 26. *Stats.txt* files, tracking statistics.

Parameter Name	Example of Value	Description
Av. Speed	44.7585	Average of <i>Instant_Speed</i>
Av. Accel	68.3577	Average of <i>Instant_Accel</i>
Mobility Rate	0.999677	Rate of <i>Instant_Speed</i> above a certain value ( <i>ana.mobs</i> )
Visible Frames	22111	Number of detections in the arena
Visible Time	884.471	(Visible Frames)/( Video Framerate)
Invisible Frames	372	(Analyzed Video Frames) - ( Visible Frames)
Invisible Time	14.91	(Visible frames)/( Video Framerate)
First Visible Frame	1	First detection frame number
Last Visible Frame	22481	Last detection frame number
Visibility Rate	0.983415	(Visible Frames)/( Visible Frames+ Visible Frames)
Invisibility Rate	0.0165853	(1)-(Visibility Rate)
Explored Areas	194	Number of areas with assignments from <i>Exploration</i>
Number of Areas	201	Size of the area matrix from <i>Exploration</i> , if normalization is selected, the image areas beyond the maximum or minimum position of the animal will not be computed.
Exploration Rate	0.967548	(Number of Areas)/(Number of Explored Areas (from <i>Exploration</i> ))
Total Distance	164834	Total Swimming distance in mm
Transitions to White	0	Number of times where the animals appears in the image (from <i>Transitions</i> )
Transitions to Black	0	Number of times when the animals disappears in the image (from <i>Transitions</i> )
Number of Frozen Events	1	Count of the Frozen events from <i>FrozenEvents</i>
Total Time Frozen	6.63131	Total time in frozen state from <i>FrozenEvents</i>
Average Time Frozen	1.11772	Average time of a frozen state from <i>FrozenEvents</i>

## Advanced parameters

### 1. Graphic parameters (*ColorIni.txt*)

These parameters have only cosmetic effects, and determine different graphical parameters of the application. When executing the software (or when “New Project” option is selected), it loads as default parameters the values located in the *ColorIni.txt* from the application folder. The graphic parameters are separated in two groups and each one has an unique code for identification. Therefore, Table 27 show general parameters related to font sizes and line width. And Table 28 defines color properties.

The width and size values are scaled in according to the image or window resolution, and the color values represent the RGB color values in inverse order (BGR), in accordance to the OpenCv representation. The graphic parameters are detailed as follows:

Table 27. *ColorIni.txt* file, general parameters.

GENERAL_PARAMETERS		
Code	Default Value	Function
traj.long	25	Length of the trajectory displayed in the tracking images.
roiL.widt	2	Line width of the rectangles showing the arena ROI.
tral.widt	1	Line width of the trajectories.
stal.widt	1	Line width of zone areas in the graphical outputs.
lbl.size	2	Label size, for the tracking images
font.size	2	Font size.

Table 28. *ColorIni.txt* file, color parameters.

COLOR_PARAMETERS		
Code	Default Value	Function
rea.bgnd	255 255 255	Color for the background of the virtual arena, for the population graphical outputs.
staL.colr	0 0 0	Color for the zone areas in the graphical outputs.
roiU.colr	255 0 0	Color for the non-selected arena ROIs
roiS.colr	0 255 0	Color for the selected arena ROI.
roNU.colr	255 0 0	Color for the non-selected arena names.
roNS.colr	0 255 0	Color for the selected arena name.
roiM.colr	0 0 255	Color for the tracking areas

## 2. Configuration parameters (*Configuration.txt*)

The configuration parameters are stored in the *Configuration.txt* files generated for each project. When executing the software (or when “New Project” option is selected), it loads as default parameters the values located in the *ConfigurationIni.txt* from the application folder. The configuration parameters are separated in groups and each one has an unique code for identification, these groups are detailed in Tables 29-59. The configuration parameters are detailed as follows:

Table 29. *Configuration.txt* files, calibration parameters.

GENERAL PARAMETERS			
Code	Default Value	Function	Special Values
exe.thre	16	Number of execution threads. Maximum number of parallel processes using for processing by the software	1/ 0 = Non-Parallel Implementation

Table 30. *Configuration.txt* files, calibration parameters.

CALIBRATION_PARAMETERS			
Code	Default Value	Function	Special Values
cal.size	20	Size of each square of the checkerboard pattern used for calibration. (Available on the Interface)	
cal.cols	10	Calibration pattern cols (internally 1 is subtracted to this value to this). (Available on the Interface)	
cal.rows	8	Calibration pattern rows (internally 1 is subtracted to this value to this). (Available on the Interface)	
cal.dist	1	Calibration distortion model. (Available on the Interface)	0=Rad3 1=Rad3+Tangent2 2=Rad3+Tangent2 3=Rad3+Tangent2+Prism4

Table 31. *Configuration.txt* files, arena definition parameters.

ARENA_DEFINITION_PARAMETERS			
Code	Default Value	Function	Special Values
roi.mode	0	Mode of arena definition selection. (Available on the Interface).	0=Automatic 1= Manual
roi.thre	150	Threshold value for selecting the tracking area (0:255). (Available on the Interface).	
roi.poly	1	The tracking area is approximated by a polygonal shape, this parameter specifies the approximation accuracy. This is the maximum distance between the original curve and its approximation. Only in automatic selection. (Available on the Interface).	
roi.elms	7	Size of the morphological element applied to in the closing operation of the tracking area. (Available on the Interface).	
roi.dilt	1	Number of times the erosion operation is applied. (Available on the Interface).	
roi.erot	4	Number of times the erosion operation is applied. (Available on the Interface).	
roi.mins	100000	Minimum area in pixels for creating an arena. Only in automatic selection. (Available on the Interface).	
roi.fite	0	Fit ellipse. Tries to convert the tracking area to the minimum enclosing circle of the selected pixels. Only for the manual selection. (Available on the Interface).	0=Disabled 1=Enabled
roi.redr	1	If previous parameter is enabled, specifies the reduction in pixels of the radio of the minimum enclosing circle (Available on the Interface).	

Table 32. *Configuration.txt* files, background parameters.

BACKGROUND_PARAMETERS			
Code	Default Value	Function	Special Values
bgs.mode	0	Mode (Available on the Interface)	0=Disabled 1=Enabled
bgs.nums	500	Sets the number of last frames that affect the background model	
bgs.thre	25	Threshold on the squared <i>Mahalanobis</i> distance between the pixel and the model to decide whether a pixel is well described by the background model. This parameter does not affect the background update.	
bgs.shad	0	Model with Shadows	0=Disabled 1=Enabled
bgs.numg	5	Number of Gaussian functions used in the model.	
bgs.ratb	0.99	Background ratio (0:1)	
bgs.lstp	1.E-06	Learning rate of the background model (0:1)	<0=Automatic 0=Not update 1=Update Always

Table 33. *Configuration.txt* files, preprocessing parameters.

PREPROCESSING_PARAMETERS			
Code	Default Value	Function	Special Values
pre.gfil	5	Gaussian filter preprocessing applied for every frame. The value is used as the size of the Gaussian filter	0=Disabled 1=Enabled
pre.norm	0	Enables or disables an image normalization preprocessing	0=Disabled 1=Enabled

Table 34. *Configuration.txt* files, detection parameters.

DETECTION_PARAMETERS			
Code	Default Value	Function	Special Values
det.type	0	Type of detection (not used and reserved for modding)	
det.thre	90	Threshold value for the segmentation (0:255). (Available on the Interface)	0=Otsu value

Table 35. *Configuration.txt* files, detection opening/closing parameters.

DETECTION_OPENING_CLOSING			
Code	Default Value	Function	Special Values
det.opcl	0	Type of morphological operation performed.	0 = Opening (Erosion+Dilation) 1 = Closing (Dilation+Erosion)
det.elms	3	The size of the structuring element for the opening/closing operation	0 = Disabled
det.dilt	2	Iterations of the dilation for the opening/closing operation	0 = Disabled
det.erot	2	Iterations of the erosion for the opening/closing operation	0 = Disabled

Table 36. *Configuration.txt* files, detection dilation/erosion parameters.

DETECTION_DILATION_EROSION			
Code	Default Value	Function	Special Values
det.erdi	0	Type of morphological operation performed.	0 = Dilation 1 = Erosion
det.elss	0	Size of the morphological element applied to in the erosion operation of the detected objects after closing.	0=Disabled
det.ertt	0	Number of times the erosion operation is applied.	0=Disabled

Table 37. *Configuration.txt* files, detection filter parameters.

DETECTION_FILTER			
Code	Default Value	Function	Special Values
det.filt	1	Use filtering on detected bodies. (Available on the Interface).	0=Disabled 1=Enabled
det.maxs	1500	Maximum size in pixels of the detected bodies (Available on the Interface).	
det.mins	150	Minimum size in pixels of the detected bodies (Available on the Interface).	
det.maxr	0	Maximum size in pixels of the minimum enclosing circle of the detected bodies. (Available on the Interface).	0=Disabled
det.minr	0	Minimum size in pixels of the minimum enclosing circle of the detected bodies. (Available on the Interface).	0=Disabled
det.mash	0	Maximum rate between the major and the minor radius of the minimum ellipse fitting the detected bodies.	0=Disabled
det.mish	0	Minimum rate between the major and the minor radius of the minimum ellipse fitting the detected bodies.	0=Disabled
det.minf	0	Minimum fill rate (ratio between the area of the minimum ellipse fitting the body and the actual number of pixels detected in the body).	0=Disabled

Table 38. *Configuration.txt* files, Kalman filter type.

KALMAN_FILTER_TYPE			
Code	Default Value	Function	Special Values
kal.mode	2	Kalman filter algorithm (not used and reserved for modding).	

Table 39. *Configuration.txt* files, Kalman filter parameters.

KALMAN_FILTER_PARAMS			
Code	Default Value	Function	Special Values
kal.time	0.25	Time increment magnitude of the Kalman filter.	
kal.pron	1.E-01	Estimated process noise magnitude (determines expected changes in acceleration) of the Kalman filter.	
kal.mean	1.E-05	Estimated measurement noise magnitude of the Kalman filter (determines how much the Kalman filter will follow the detected positions).	
kal.errc	1.E-01	Magnitude to initialize the a posteriori error covariance matrix of the Kalman filter.	

Table 40. *Configuration.txt* files, Kalman filter acceptance parameters.

KALMAN_FILTER_ACCEPTANCE			
Code	Default Value	Function	Special Values
kal.disf	50	Frame distance condition, defines the maximum allowed pixels the animal can displace between consecutive frames. (Available on the Interface).	
kal.sich	0.4	Size change condition.	kal.sich

Table 41. *Configuration.txt* files, Kalman filter first delete condition parameters.

KALMAN_FILTER_DELETE1			
Code	Default Value	Function	Special Values
kal.dund	1	Maximum number of consecutive frames, a track has not been assigned to any detection, before mark this track as inactive.	

Table 42. *Configuration.txt* files, Kalman filter second delete condition parameters.

KALMAN_FILTER_DELETE2			
Code	Default Value	Function	Special Values
kal.dage	10	Minimum size of a track to be used by the algorithm, tracks smaller than this size will be deleted.	
kal.dmax	8	Minimum number of detections for a track to be considered valid.	

Table 43. *Configuration.txt* files, Kalman filter number of tracks.

KALMAN_FILTER_NUMBEROFTRACKS			
Code	Default Value	Function	Special Values
kal.ntra	1	Defines the number of animals per arena. (Available on the Interface).	

Table 44. *Configuration.txt* files, identity algorithm.

KALMAN_MULTITRACKING_IDENTITY_ALGORITHM			
Code	Default Value	Function	Special Values
kal.idal	0	Selects track identification algorithm. (Available on the Interface).	0=Not Id 1=Hist Only 2=Hist + TCM
kal.fdis	0	Frame distance condition, defines the maximum allowed pixels the animal can displace between consecutive frames, to consider two tracks compatibles.	0=Calculated automatically

Table 45. *Configuration.txt* files, identity algorithm, comparison parameters.

KALMAN_MULTITRACKING_COMPARISON			
Code	Default Value	Function	Special Values
kal.cmse	0.2	Minimum size change between to samples of a track to compute the identity matrix	
kal.corH	1	Use Histogram Correlation Distributions	0=Disabled 1=Enabled
kal.cmhc	0.7	Minimum histogram correlation between to samples of a track to compute the identity matrix (For CMT algorithm)	
kal.corN	100	Number of classes in the Correlation Distributions	
kal.shaC	10	Maximum Shape alignment error (For Shape algorithms)	
kal.shaN	500	Number of classes in the Shape Distributions (used for optimization only)	

Table 46. *Configuration.txt* files, identity algorithm, evaluation parameters.

KALMAN_MULTITRACKING_EVALUATION			
Code	Default Value	Function	Special Values
kal.mcnp	500	Maximum number of comparisons, tries to limit the number of sample comparisons used for two tracks	
kal.mAvg	1	Use Avg. Use a mean metric instead of a max similarity value	0=Disabled 1=Enabled
kal.gStd	0.05	Std. Dev of the Gaussian function ruling the decay of the weight in the Histogram Correlation Distribution	

Table 47. *Configuration.txt* files, identity algorithm, selection parameters.

KALMAN_MULTITRACKING_SELECTION			
Code	Default Value	Function	Special Values
kal.rGrp	2	Restriction to groups, Allow the first group/all groups/all tracks to be matched	0=All tracks 1>All groups 2=First group
kal.hOrd	1	Use High order Correlation. Takes in account the similarity with all fragments, without reducing the speed.	

Table 48. *Configuration.txt* files, identity algorithm, feature parameters.

KALMAN_MULTITRACKING_FEATURE_PARAMETERS			
Code	Default Value	Function	Special Values
kal.hiss	20	Size (number of color clusters) used by the histograms. (Available on the Interface).	
kal.tcmr	25	Max distance from the body center used by the TCM maps. (Available on the Interface).	
kal.tcnd	0	Selects the data type the TCM maps will use. (Available on the Interface).	0=8 Bits 1=16 Bits 2=32 Bits
kal.hist	500	History, maximum number of samples used by every track in the algorithm. (Available on the Interface).	

Table 49. *Configuration.txt* files, tracking collision parameters.

KALMAN_MULTITRACKING_COLLISION_PARAMETERS			
Code	Default Value	Function	Special Values
kal.advr	0.8	Distance modifier, will define the minimum advantage for two tracks with the same closest detection as $kal.disf \times kal.advr$ for detecting a collision.	
kal.advm	10	Minimum advantage value in pixels for detecting a collision between two tracks.	
kal.cnft	20	Number of frames conflicted tracks is kept in memory. A track is marked as conflicted if it collides with another track. The fusion algorithm will use these tracks to try to solve some collisions.	

Table 50. *Configuration.txt* files, tracking fusion parameters.

KALMAN_MULTITRACKING_FUSSION_PARAMETERS			
Code	Default Value	Function	Special Values
kal.tfmi	5	Minimum track age to be fused.	
kal.tfma	10	Maximum track age to be fused.	
kal.tdma	10	Maximum temporal distance between the active track and the suitable track for fusion.	
kal.acor	0.6	Minimum mean histogram correlation value between the active track and the suitable track for fusion.	
kal.bcor	0.5	Minimum best histogram correlation value between the active track and the suitable track for fusion.	

Table 51. *Configuration.txt* files, minimum track size.

KALMAN_MULTITRACKING_TRACK_PARAMETERS			
Code	Default Value	Function	Special Values
kal.mins	50	Minimum track size for a track to be marked as a long track, used in the track identification algorithm.	

Table 52. *Configuration.txt* files, identity algorithm correlation parameters.

KALMAN_MULTITRACKING_CORRELATION_PARAMETERS			
Code	Default Value	Function	Special Values
kal.idgb	0	Minimum correlation value to be accepted by the Hungarian algorithm in the track identification algorithm. The algorithm will still be able to identify the tracks if only one of the assignments has a smaller correlation value. (Available on the Interface).	
kal.idla	0	Minimum average correlation value to be accepted to identify one single long track in the track identification algorithm. (Available on the Interface).	
kal.idlb	0	Minimum best correlation value to be accepted to identify one single long track in the track identification algorithm. (Available on the Interface).	
kal.idsa	0	Minimum average correlation value to be accepted to identify one single short track in the track identification algorithm. (Available on the Interface).	
kal.idsb	0	Minimum best correlation value to be accepted to identify one single short track in the track identification algorithm. (Available on the Interface).	

Table 53. *Configuration.txt* files, identity algorithm, maximum number of tracks.

KALMAN_MULTITRACKING_OTHER_PARAMETERS			
Code	Default Value	Function	Special Values
kal.idff	500	Maximum number of tracks, kept in memory for the track identification algorithm. When the number of tracks reaches this number, the algorithm will try to identify the tracks to release memory, before continuing the tracking.	

Table 54. *Configuration.txt* files, output parameters.

OUTPUT_PARAMETERS			
Code	Default Value	Function	Special Values
out.step	10	Allows changing the number of frames in the tracking processing stage to update the interface or reaching a safe processing point (check for stop signal).	
out.wind	0	Enables to show the tracking results in real time, showing a free resizable window of the current tracking for every arena. This option is useful for demonstration or debugging, but slows down significantly the processing speed.	0=None 1=Tracking 2=Trajectory
out.ftxt	1	Modes of txt output. Changes the amount of plain text files the software generates in the results.	0=None 1=Main stats 2>All text files
out.fjpg	1	Enables to save an image file for every arena and frame in real time showing the tracking results. This option is useful for demonstration or debugging, but slows down significantly the processing speed. Also allows enabling or disabling the spatial stats image results.	0=None 1=Spatial stats only 2=Tracking 3=Trajectory 4=Trajectory
out.pnam	TestProject	Project name, it will be used as a part of the filenames of the output files.	

Table 55. *Configuration.txt* files, data analysis arena parameters.

DATA_ANALYSIS_arena			
Code	Default Value	Function	Special Values
ana.norm	0	Normalize arenas. If enabled, the zones will be computed according to the extremes of the detected positions instead of using the full image.	0=Disabled 1=Enabled
ana.aror	3	If the arenas are not symmetrical, and do not have the same orientation. This button allows changing the orientation of the arenas, in the projected virtual arena for the entire population.	0=Same Orientation 1=Horizontal Mirror 2=Vertical Mirror 3=Vertical & Horizontal Mirror

Table 56. *Configuration.txt* files, data analysis zone parameters.

DATA_ANALYSIS_ZONE			
Code	Default Value	Function	Special Values
ana.nzon	30	Max number of areas in the edge spatial results.	
ana.zsizq	50	Distance in mm of each area in the spatial results.	

Table 57. *Configuration.txt* files, data analysis speed parameters.

DATA_ANALYSIS_SPEED			
Code	Default Value	Function	Special Values
ana.spsa	2	Sampling distance (in frames) to estimate the instantaneous speeds.	
ana.mobs	1	Speed threshold to estimate the mobility rate of the animal in the stats.	

Table 58. *Configuration.txt* files, data analysis frozen events parameters.

DATA_ANALYSIS_FROZEN			
Code	Default Value	Function	Special Values
ana.fmmf	5	A frozen event is detected when the animal has moved less than this value (in mm) in a certain amount of time (ana.ftim).	
ana.ftim	3	A frozen event is detected when the animal has moved less than ana.fmmf in this amount of time (in seconds).	

Table 59. *Configuration.txt* files, data analysis transitions parameters.

DATA_ANALYSIS_TRANSITIONS			
Code	Default Value	Function	Special Values
ana.ttim	7	A transition is detected when the time between two consecutive detections exceed this value (in seconds).	

Table 60. *Configuration.txt* files, data analysis post process parameters.

DATA_ANALYSIS_POSTPROCESS			
Code	Default Value	Function	Special Values
ana.inte	0	Interpolate holes in the trajectory, using a linear interpolation algorithm. This parameter can be changed without redoing the analysis.	0=Disabled 1=Enabled
ana.intf	25	Maximum size of a trajectory hole (in frames) where interpolation will be applied. This parameter can be changed without redoing the analysis.	
ana.smoo	0	Use a moving average to smooth trajectories. This parameter can be changed without redoing the analysis.	0=Disabled 1=Enabled

Table 61. *Configuration.txt* files, data analysis other parameters.

DATA_ANALYSIS_OTHER			
Code	Default Value	Function	Special Values
ana.rvis	0.05	If visibility rate is smaller than this value, normalization will not be applied even if the option is selected.	

Table 62. *Configuration.txt* files, main video parameters.

MAIN_VIDEO_PARAMETERS			
Code	Default Value	Function	Special Values
oth.mini	0	Starting point in minutes (for all sequences) for analysis.	
oth.mend	1.E+10	Ending point in minutes (for all sequences) for analysis.	
oth.atyp	8	Arena type (not used and reserved for modding).	
oth.frat	25	Frame Rate (not used and reserved for modding).	
oth.rees	1	Change of image size (not used and reserved for modding).	

## References

1. Zhang, Z. A flexible new technique for camera calibration. *IEEE Trans. Pattern Anal. Mach. Intell.* **22**, 1330–1334 (2000).
2. Bouguet, J.-Y. MATLAB calibration tool. (2014).
3. Zivkovic, Z. Improved adaptive Gaussian mixture model for background subtraction. in *Pattern Recognition, 2004. ICPR 2004. Proceedings of the 17th International Conference on* **2**, 28–31 (IEEE, 2004).
4. Zivkovic, Z. & van der Heijden, F. Efficient adaptive density estimation per image pixel for the task of background subtraction. *Pattern Recognit. Lett.* **27**, 773–780 (2006).
5. Welch, G. & Bishop, G. *An Introduction to the Kalman Filter*. (Department of Computer Science, University of North Carolina , 2006).
6. Kuhn, H. W. The Hungarian method for the assignment problem. *Nav. Res. Logist. Q.* **2**, 83–97 (1955).